

## 1. BF-48AS Specification

	B F - 4 8 A S
Injection volume	4 $\mu$ l or 2 $\mu$ l
Sample volium	30 $\mu$ l to 500 $\mu$ l (variable)
Number of samples	48
Dilution	Undiluted, or by 2 to 100 times
Dilute solution	Rinse solution
Power source	AC100V, 50/ 60Hz, max3A
Size	200(W) $\times$ 420(D) $\times$ 485(H)mm
Weight	14.3kg
Hygrothermal environment	15 to 35, 20 to 80%RH, no condensation